

# CURVEMATE



# CURVEMATE

FOR A HAPPY MIND  
AND HEALTHY LIFE

## Caution :

- Avoid exerting strong pressure on the head.
- Avoid use in high humidity.
- Keep in dry places.

## I N S T R U C T I O N


### How to set a battery :

Operates on two AA (UM-3) batteries

If used continuously :

Regular battery lasts/approx. 90 minutes

Alkaline battery lasts/approx. 180 minutes




Set ⊕ ⊖ of the battery on  
right order.


Exert light pressure on  
this part of open and  
close the battery chamber  
cover.

### CURVEMATE :

Curvemate is a handy machine designed to  
ease muscular tension. With Curvemate, you  
can relieve your stiff shoulders and tired legs  
instantly and conveniently. So keep  
it handy wherever  
you go!



**Instructions :** To operate curvemate, simply press the head lightly against desired parts of your body.



Put cap back  
On after use.



For relief of neck and shoulder  
stiffness



For shaping up



For relief of tired muscles