



**BATTERY OPERATED
WITH ON/OFF SWITCH**

- Acts as a lumbar support for the lower back
- Soothes tired backs and stiff necks
- Gentle massage action
- Helps relieve the lower back pains that occur during the later months of pregnancy*
- Easy to use anywhere
- Soft corduroy shell with thick Polyurethane foam filler

CHIROPRACTIC MASSAGE PILLOW



BATTERY INSTRUCTIONS

1. Open flap on back of pillow.
2. Do not remove entire vibrator unit.
3. Insert two D-size batteries in vibrator unit.
4. After batteries are installed close pillow & enjoy.
5. Push red button at end to start and stop vibration

Batteries not included/takes 2 D batteries

